

DSM Biomedical CarboSil® 20 90A Thermoplastic Silicone Polycarbonate Polyurethane (TSPCU)

Category : Polymer , Thermoplastic , Polycarbonate (PC) , Polyurethane, TP , Silicone Polyurethane, Polycarbonate Based

Material Notes:

A ground-breaking copolymer that combines the biocompatibility and biostability of conventional silicone elastomers with the processability and toughness of thermoplastic polycarbonate-urethanes. The silicone portion of CarboSil® TSPCU works synergistically with the polycarbonate component to improve stability. This medical grade polymer is highly biocompatible and well suited to be used in many types of medical devices. Flexible CarboSil® TSPCU is comparable in tensile strength to traditional polycarbonate urethanes and derives additional biostability from the silicone portion. The material is adaptable to various fabrication techniques to accommodate many different device shapes. It can be extruded, injection or compression molded, solvent bonded, dipped coated and sprayed. Widely Used Because of its tensile strength, biocompatibility and flexibility, CarboSil® TSPCU is used in a wide range of medical applications, including cardiovascular and nervous system electrostimulation, continuous glucose monitoring, drug eluting and orthopedic implants. Tailor Made CarboSil® TSPCU can be enhanced with SME® technology to incorporate end groups that can address the needs of specific device applications (silicone end groups are standard). This eliminates the need for additional surface processing steps after the device component is fabricated. Summary of Product Benefits Biostable and biocompatible Adaptable to many different processing techniques Excellent mechanical properties Thromboresistant Enhanced with SME® technology Comprehensive FDA Master File Grade denotes silicone content (%) Information provided by DSM Biomedical.

Order this product through the following link:

http://www.lookpolymers.com/polymer_DSM-Biomedical-CarboSil-20-90A-Thermoplastic-Silicone-Polycarbonate-Polyurethane-TSPCU.php

Physical Properties	Metric	English	Comments
Specific Gravity	1.16 g/cc	1.16 g/cc	ASTM D792
Melt Flow	27 g/10 min @Load 1.20 kg, Temperature 224 °C	27 g/10 min @Load 2.65 lb, Temperature 435 °F	ASTM D1238

Mechanical Properties	Metric	English	Comments
Hardness, Shore A	90	90	ASTM D2240
Tensile Strength, Ultimate	39.98 MPa	5798 psi	ASTM D1708
Tensile Strength, Yield	7.460 MPa @Strain 50.0 %	1082 psi @Strain 50.0 %	ASTM D1708
	9.832 MPa @Strain 100 %	1426 psi @Strain 100 %	ASTM D1708
	26.06 MPa @Strain 300 %	3779 psi @Strain 300 %	ASTM D1708

Elongation at Break Mechanical Properties	424 % Metric	424 % English	ASTM D1708 Comments
Flexural Modulus	0.0407 GPa	5.91 ksi	ASTM D790
Tear Strength	87.7 kN/m	500 pli	Die 'C'; ASTM D624
Taber Abrasion, mg/1000 Cycles	57	57	H-18 Wheel; ASTM D1044
Compression Set	15 % @Temperature 25.0 °C	15 % @Temperature 77.0 °F	22 hours; ASTM D395

Thermal Properties	Metric	English	Comments
Glass Transition Temp, Tg	-4.00 °C	24.8 °F	DSC

Processing Properties	Metric	English	Comments
Processing Temperature	190 - 210 °C	374 - 410 °F	Optimum Extrusion Conditions

Descriptive Properties	Value	Comments
26 week Carcinogenicity Study in the Transgenic ras H2 Mouse Model	No increase in induced tumor formation	
Chronic Toxicity, Subcutaneous Implant	No evidence of systemic toxicity	
Color	Clear to Amber Pellets	
Cytotoxicity	No evidence of causing cell lysis or toxicity	
Genotoxicity: Bacterial Reverse Mutation (95% ethanol extract)	Non-mutagenic	
Genotoxicity: Bacterial Reverse Mutation (saline extract)	Non-mutagenic	
Genotoxicity: In Vitro Chromosomal Aberration	Non-genotoxic	
Hemolysis	Non-hemolytic	
ISO Intracutaneous Irritation	No evidence of significant irritation	
ISO Maximization Sensitization	No evidence of causing delayed dermal contact sensitization	
Mouse Bone Marrow Micronucleus	Non-genotoxic	
Muscle Implantation, 12 weeks	Non-irritant	
Muscle Implantation, 2 weeks	Non-irritant	
USP and ISO Systemic Toxicity	No evidence of systemic toxicity	
USP Pyrogen Study	Non-pyrogenic	

Contact Songhan Plastic Technology Co.,Ltd.

Website : www.lookpolymers.com

Email : sales@lookpolymers.com

Tel : +86 021-51131842

Mobile : +86 13061808058

Skype : lookpolymers

Address : United North Road 215,Fengxian District, Shanghai City,China